

To whom it may concern

Shelf-life and storage of Nidacon products.

The following products have been extensively tested for both shelf-life and storage conditions. The same shelf-life applies even after opening for all products if aseptic procedures are being used.

Product	Shelf-life	Storage unopened	Storage opened
PureSperm100	24 months	2-27°C	2-8°C
PureSperm 40/80/90	24 months	2-27°C	2-8°C
PureSpermBuffer	24 months	2-40°C	2-8°C
ProInsert	24 months	17-27 °C	17-27 °C
SpermVitalStain	24 months	10-40°C	17-27 °C
Nidoil	24 months	2-27°C	2-8°C
PureSpermWash	12 months	2-30°C	2-8°C
SpermCatch	12 months	2-30°C	2-8°C
SpermCryoProtec	12 months	2-30°C	2-8°C
PureSpermSpeediKit	12months	2-30°C	2-8°C
VitriBlast	12 months	2-30°C	2-8°C
ThermoBlast	12months	2-30°C	2-8°C

All documents and records of test data have been reviewed by the relevant authorized persons.

Signature:



Marina Danilova
QA Coordinator
Nidacon International AB